

ENT

Functional Endoscopic Sinus Surgery (FESS)



I N K A™

SURGICAL INSTRUMENTS

DISTRIBUTED BY:



SURGICAL ASIA Sdn Bhd (1272656-H)

Suppliers of Fine Quality Instruments

No. 8, Lebuhr Perusahaan Klebang 9A,
Kawasan Perusahaan IGB, 31200 Chemor, Perak, Malaysia

Phone: +605 210 8080 • www.inkasurgical.asia



Exclusive Distributor

OBEX Medical Ltd

0800 656 239

www.inkasurgical.co.nz



we thank you for your support

Disclaimer: Product availability is subject to regulatory approval in individual markets.
Not all products are available in all markets. Contact us to find out more.

BRC-00077


contents

01	Legend	42	Antrum - Frontal Sinus Punches
03	Nasal Packing Forceps	44	Trocar - Cannulas
04	Specula - Nasal	45	Probes
05	Specula - Septum	46	Knives
06	Nasal Scissors	48	Curettes
07	Delicate Scissors - Forceps	51	Curettes - Dissectors
10	Forceps	52	Elevators
15	Through Cutting Forceps	54	Cannulas - Irrigators
20	Cup Forceps		
22	Giraffe Forceps		
28	Grasping - Cutting Forceps		
29	Biopsy Forceps		
32	Backbiter Punches		
36	Rotating Backbiter		
37	Terrier Backbiter		
38	Sidebiting Punches		
39	Sphenoidal Punches		

legend

Reverse Action Handles


Reverse Action Handles are thumb activated and reduces distal tip movement.



Take-apart joint



Giraffe Forceps



Vertical action
in line with shaft



Horizontal action

Connectors




SUTTER
Standard connector for all
Bipolar Forceps

- A EU flat connector
- B US2-pin connector


legend

SuperGliss® Tips




Extra Low Profile -

ELP tips are shorter and more delicate. Ideal for fine microsurgical interventions.




Micro Tip -

Makes insulation disappear from the surgeon's sight and opens up the view through the tips.




Classic Tip -

Plateau shape for strong grip. Allows dissection, grasping and coagulation of larger structures.





TEO Tip -

Ideal for endoscope assisted procedures. Short round tips provide safety and better grip.



Zhora Tip -

Very delicate tips. Can be clearly seen in the operating room due to the 10° upward skew. Narrow and confined spaces are advantageous.




SuperGrip -

Pinpoint grasping and better grip due to microstructured tips and clearly defined edges.


Colour	Tip Size	EU Flat Connector	US2-pin Connector
Yellow	0.2mm		
Blue	0.4mm		
Red	0.7mm		
Green	1mm		
Black	1.2mm+		

nasal packing forceps

Instruments @ 50% · Inserts to Scale



LISTER
Sinus Forceps
straight
12.5cm
13620.12



TILLEY
Nasal Dressing Forceps
17.5cm
13624.17



TILLEY
Nasal Dressing Forceps
slender pattern
14cm
13625.14

specula - nasal

Instruments @ 50% · Inserts to Scale



TIECK-HALLE
Nasal Specula
for infants
13.5cm

A 33112.01 Small
B 33112.02 Large




VIENNA
Nasal Specula
14cm

D 33116.01 Small
E 33116.02 Medium
F 33116.03 Large


specula - septum

Instruments @ 50% · Inserts to Scale




nasal scissors

Instruments @ 50% · Inserts to Scale



Nasal Scissors
with tubular shaft,
serrated blades,
11mm cutting edge

- A **33195.01** Straight
- B **33195.03** Left curved
- C **33195.04** Right curved





Nasal Scissors
with 11mm smooth blades,
sliding shaft

- A **33197.01** Straight
- B **33197.03** Left curved
- C **33197.04** Right curved



delicate scissors - forceps

Instruments @ 50% · Inserts to Scale





STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm,
long shaft

- A 33257.25/13** Straight
- B 33257.60/13** 60° up-curved
- C 33257.90/13** 90° up-turned



STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm,
x-long shaft

- 33257.25/16** Straight
- 33257.60/16** 60° up-curved
- 33257.90/16** 90° up-turned




BELLUCCI
Paediatric Nasal Scissors
delicate pattern,
straight blades 5mm

- 31650.13/13** Shaft 13cm
- 31650.13/15** Shaft 15cm


delicate scissors - forceps

Instruments @ 50% · Inserts to Scale




HARTMANN
Alligator Forceps
serrated

30420.14 Shaft 14cm
30420.16 Shaft 16cm




HARTMANN
Round Cup Forceps
cup ø 2.5mm, straight

30445.25-12 Long shaft, 12.5cm
30445.25-15 X-long shaft, 15cm





WULLSTEIN
Micro Round Cup Forceps
cup ø 1mm, straight

31588.10/13 Long shaft, 13cm
31588.10/15 X-long shaft, 15cm



delicate scissors - forceps

Instruments @ 50% · Inserts to Scale





STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm,
xx-long shaft

33257.25/19 Straight
33257.60/19 60° up-curved
33257.90/19 90° up-curved


forceps

Instruments @ 50% · Inserts to Scale





forceps

Instruments @ 50% · Inserts to Scale



WIGAND
Nasal Forceps
for grasping polyps,
60° up-curved, round tip
cup jaws 14x6mm

33269.60





TILLEY HENCKEL
Ethmoid Forceps
with graduated shaft,
straight

- A 33270.45** 4.5mm bite
- B 33270.55** 5.5mm bite
- C 33270.65** 6.5mm bite
- D 33270.80** 8mm bite

forceps

Instruments @ 50% · Inserts to Scale





TAKAHASHI
Ethmoid Forceps

Straight

- A 21617.02** 2.5mm bite
- B 21617.04** 4mm bite



Up-curved

- A 21618.02** 2.5mm bite
- B 21618.04** 4mm bite



BLAKESLEY
Nasal Suction Forceps
with cleaning stilette,
straight,

- A 33324.25** 2.5mm bite, Ø 4mm
- B 33324.35** 3.5mm bite, Ø 4mm
- C 33324.40** 4mm bite, Ø 5mm




BLAKESLEY
Nasal Suction Forceps
with cleaning stilette,
30° up-curved,

- A 33325.25** 2.5mm bite, Ø 4mm
- B 33325.35** 3.5mm bite, Ø 4mm
- C 33325.40** 4mm bite, Ø 5mm


forceps

Instruments @ 50% · Inserts to Scale



WEIL-BLAKESLEY
Nasal Suction Forceps
with Luer-Lock connector

33251.48 4.8mm bite, straight
33254.48 4.8mm bite, up-curved



TAKAHASHI
Nasal Suction Forceps
with cleaning stilette

Straight



A 33328.25 2.5mm bite
B 33328.35 3.5mm bite

30° Up-curved

A 33329.25 2.5mm bite
B 33329.35 3.5mm bite




forceps

Instruments @ 50% · Inserts to Scale



BLAKESLEY
Frontal Suction Forceps
with Luer-Lock connector,
cleaning stilette, straight,
2.5mm bite

33334.25



BLAKESLEY
Frontal Suction Forceps
with Luer-Lock connector,
45° up-curved, 2.5mm bite,
cleaning stilette

33335.25

through cutting forceps


Instruments @ 50% · Inserts to Scale



through cutting forceps

Instruments @ 50% · Inserts to Scale

			Outside Width	Inside Cut
		A	3.3mm	2mm
		B	3.9mm	2.5mm
		C	4.6mm	3mm
		D	5.1mm	3.5mm
		E	5.7mm	4mm





WEIL-BLAKESLEY
Through cutting Nasal Forceps
with serrated lower jaw,
straight

- A 33261.20** 2mm cut
- B 33261.25** 2.5mm cut
- C 33261.30** 3mm cut
- D 33261.35** 3.5mm cut
- E 33261.40** 4mm cut


through cutting forceps

Instruments @ 50% · Inserts to Scale



through cutting forceps

Instruments @ 50% · Inserts to Scale




WEIL-BLAKESLEY
Through Cutting Nasal Forceps
with suction

Straight

- A 33263.25** 2.5mm bite
- B 33263.36** 3.6mm bite

Up-curved

- A 33266.25** 2.5mm bite
- B 33266.36** 3.6mm bite



GRUENWALD
Nasal Cutting Forceps with
Luer-Lock cleaning channel,
through cutting

Straight



- C 33338.20** 2mm cut
- D 33338.25** 2.5mm cut
- E 33338.30** 3mm cut

45° Up-curved

- C 33339.20** 2mm cut
- D 33339.25** 2.5mm cut
- E 33339.30** 3mm cut



through cutting forceps

Instruments @ 50% · Inserts to Scale





cup forceps

Instruments @ 50% · Inserts to Scale





Sinus Cup Forceps
round cups, $\varnothing$ 3.5mm
33241.35



Round Cup Forceps
tubular shaft, cup $\varnothing$ 4mm,
straight
33928.16 16cm
33928.18 18cm



cup forceps

Instruments @ 50% · Inserts to Scale



Biopsy & Grasping Forceps
oval cups 3x5mm,
tubular shaft

33924.12 12.5cm
33924.15 15cm





Biopsy & Grasping Forceps
tubular shaft, oval cups 3mm

33925.15 15cm
33925.18 18cm


giraffe forceps

Instruments @ 50% · Inserts to Scale





giraffe forceps

Instruments @ 50% · Inserts to Scale



Short Biopsy & Grasping Forceps, double-action cupped jaws 3x6mm, Luer-Lock cleaning channel, straight

33941.03





Short Biopsy & Grasping Forceps double-action cupped jaws, 3x6mm, vertical opening, Luer-Lock cleaning channel

33941.07 70° up

33941.11 110° up



giraffe forceps

Instruments @ 50% · Inserts to Scale





giraffe forceps

Instruments @ 50% · Inserts to Scale




giraffe forceps

Instruments @ 50% · Inserts to Scale






giraffe forceps

Instruments @ 50% · Inserts to Scale





grasping - cutting forceps

Instruments @ 50% · Inserts to Scale








biopsy forceps

Instruments @ 50% · Inserts to Scale




biopsy forceps

Instruments @ 50% · Inserts to Scale




biopsy forceps

Instruments @ 50% · Inserts to Scale




70°




SKINNER
Maxillary Biopsy Forceps
double-action cupped jaws
3mm, 70° down-curved,
Luer-Lock cleaning channel,
vertical opening

33984.07



110°




SKINNER
Maxillary Biopsy Forceps
double-action cupped jaws,
3mm, 110° down-curved,
Luer-Lock cleaning channel,
horizontal opening

33984.11



backbiter punches

Instruments @ 50% · Inserts to Scale





OSTRUM
Backbiter Antrum Punches
paediatric size, 1.5mm cut

- A 33345.15** Straight
- B 33347.15** Left cutting
- C 33348.15** Right cutting




OSTRUM
Backbiter Antrum Punches
standard size, 2mm cut

- A 33345.20** Straight
- B 33347.20** Left cutting
- C 33348.20** Right cutting





OSTRUM
Backbiter Antrum Punches
large size, 2.5mm cut

- A 33345.25** Straight
- B 33347.25** Left cutting
- C 33348.25** Right cutting




backbiter punches

Instruments @ 50% · Inserts to Scale





OSTRUM
Backbiter Antrum Punches
upward cutting

- A 33350.15** 1.5mm cut
- B 33350.20** 2mm cut
- C 33350.25** 2.5mm cut




OSTRUM-STAMMBERGER
Backbiter Antrum Punch,
backward cutting, 2.5mm cut
33351.25




OSTRUM
Backbiter Antrum Punches
with backward cutting,
Luer-Lock cleaning channel,
paediatric size, cut 1.5 x 5mm

- A 33359.15** Straight
- B 33360.15** Left cutting
- C 33361.15** Right cutting




backbiter punches

Instruments @ 50% · Inserts to Scale



OSTRUM
Backbiter Antrum Punches
with backward cutting,
Luer-Lock cleaning channel,
cut 2.5 x 7mm

- A 33359.25** Straight
- B 33360.25** Left cutting
- C 33361.25** Right cutting




OSTRUM
Antrum Suction Punches
backward cutting

- A 33364.20** Cut 2mm, shaft $\varnothing$ 4mm
- B 33364.25** Cut 2.5mm, shaft $\varnothing$ 5mm



rotating backbiter

Instruments @ 50% · Inserts to Scale





OSTRUM
Retro Antrum Punches
360° rotatable, 1.5mm cut,
paediatric size

33352.15




OSTRUM
Retro Antrum Punches
360° rotatable, 2mm cut,
standard size

33352.20




BACKBITER
Rotating Antrum Punches
360° rotatable, 1.5 x 5mm cut,
Luer-Lock cleaning channel,
paediatric size

33354.15





rotating backbiter

Instruments @ 50% · Inserts to Scale




BACKBITER
Rotating Antrum Punches
360° rotatable, 2.5 x 7mm cut,
Luer-Lock cleaning channel,
standard size

33354.25




BACKBITER
Rotating Antrum Punches
Luer-Lock cleaning channel,
small cut 2 x 3.5mm

A 33355.20 Shaft up-curved
B 33356.20 Shaft down-curved





terrier backbiter

Instruments @ 50% · Inserts to Scale




OSTRUM
Terrier Antrum Punches
backward cutting, paediatric
size, small oval cut 2mm

- A 33342.20** Down-curved
- B 33343.20** Up-curved




OSTRUM
Terrier Antrum Punches
backward cutting,
large oval cut 3.8mm

- A 33342.38** Down-curved
- B 33343.38** Up-curved





sidebiting punches

Instruments @ 50% · Inserts to Scale




Sidebiting Antrum Punches through cutting with Luer-Lock cleaning channel, standard size, 3.5mm cut

- A 33372.35** Left cutting
- B 33373.35** Right cutting




Sidebiting Antrum Punches through cutting with Luer-Lock cleaning channel, paediatric size, 2.5mm cut

- A 33372.25** Left cutting
- B 33373.25** Right cutting




Sidebiting Antrum Punches through cutting with Luer-Lock cleaning channel, 40° down-biting, 3.5mm cut

- 33372.35-40** Left cutting
- 33373.35-40** Right cutting




sphenoidal punches

Instruments @ 50% · Inserts to Scale




HAJEK-KOFLER
Sphenoidal Punch
fixation nut, 360° rotatable,
3mm bite

33397.03




MC KENTY
Sphenoidal Punch
through cutting, fixation nut,
360° rotatable, 2mm cut

33400.02





sphenoidal punches

Instruments @ 50% · Inserts to Scale




HAJEK
Sphenoidal Punch
40° forward cutting,
4mm bite

33403.04




HAJEK
Sphenoidal Punch
40° backward cutting,
4mm bite

A 33406.04 Up-biting
B 33407.04 Down-biting




sphenoidal punches

Instruments @ 50% · Inserts to Scale





antrum - frontal sinus punches

Instruments @ 50% · Inserts to Scale




Antrum Punch
with stilette, bone collector,
through cutting, 2mm cut

33414.02




Antrum Punch
with stilette, bone collector,
through cutting, 3mm cut

33415.03



antrum - frontal sinus punches

Instruments @ 50% · Inserts to Scale



Frontal Sinus Punches
with articulated shaft, 70°
up-turned, 40° forward
cutting with lock,
14cm working length

- A 33417.02** 2mm bite
- B 33417.03** 3mm bite




Spare Stilette only
for removal of bone material,
20cm

33414.99


trocar - cannulas

Instruments @ 50% · Inserts to Scale



Sinus Trocar & Cannula
bill-shaped tip, O.D. 5mm,
for telescope $\varnothing$ 4mm, 16cm

33854.05



Sinus Trocar & Cannula
oblique tip, for telescope,
O.D. 5mm, $\varnothing$ 4mm, 16cm

33850.05



Sinus Trocar & Cannula
oblique tip, for telescope,
 $\varnothing$ 2.7mm, O.D. 3.3mm, 16cm

33851.33

probes

Instruments @ 50% · Inserts to Scale



Frontal Ostium Seeker Probe
double-ended, one tip 90°
bent, 21cm

33866.21



Maxillary Ostium Seeker Probes
double ended, regular pattern,
ball tips 1.2 / 2mm, 21cm

33870.21




Maxillary Ostium Seeker Probes
double ended, large pattern,
ball tips 2 / 2.5mm, 21cm

33870.22


knives

Instruments @ 50% · Inserts to Scale



Antrum Sickle Knives
19cm

- A 33875.19** Pointed tip
- B 33876.19** Blunt tip



Antrum Sickle Knives
30° angled blade, pointed tip,
19cm

- L 33875.30** Left angled
- R 33875.31** Right angled



Paediatric Sickle Knife
octagonal handle
16cm

31271.16



DRAF
Sickle Knives
pointed tip
19cm

- 33901.19** Straight shaft
- 33906.19** Bayonet shaft


knives

Instruments @ 50% · Inserts to Scale




DRAF
Flap Knives
oval blade
19cm

33902.19 Straight shaft
33907.19 Bayonet shaft



DRAF
Round Knives
45° angled, straight shaft,
19cm

A 33903.15 Ø 1.5mm
B 33903.20 Ø 2mm



DRAF
Round Knives
45° angled, bayonet shaft,
19cm

A 33908.15 Ø 1.5mm
B 33908.20 Ø 2mm




Hook Knife
Retrograde, 4mm,
16.5cm

G41-1175


curettes

Instruments @ 50% · Inserts to Scale




Antrum Curettes
oval cups, 19cm

- A 33880.01** Size 1, 7mm
- B 33880.02** Size 2, 8mm
- C 33880.03** Size 3, 9mm




Antrum Curettes
straight shaft, oblong cup,
19cm

- A 33881.02** 2mm
- B 33881.03** 2.5mm



Antrum Curettes
oblong cup, angled 30°,
19.5cm

- A 33886.02** Small, 2mm
- B 33886.03** Large, 3mm



Antrum Curettes
oblong cup, x-large,
angled 45°, 4mm,
19.5cm

33887.19

curettes

Instruments @ 50% · Inserts to Scale



5mm

Antrum Curettes
round cup, angled 45°,
19.5cm

33884.19



A

B

2mm

Frontal Sinus Curettes
oblong cup, small size,
19cm

A 33890.45 45° curved

B 33890.90 90° curved



A

B

3mm

Frontal Sinus Curettes
oblong cup, regular size,
19cm

A 33891.55 55° curved

B 33891.90 90° curved



Ø 2.8mm
sharp,
90° tips


Ø 2.3mm
semi sharp,
45° tips

Antrum Curettes
double-ended, oval cups,
90°/45° angled tips
19cm

33899.19


curettes

Instruments @ 50% · Inserts to Scale




Antrum Curettes
forward up-cutting,
small size,
19cm

33895.18




Antrum Curettes
backward down-cutting
small size,
19cm

33896.18



Antrum Curettes
forward up-cutting,
regular size,
19cm

33895.19





Antrum Curettes
backward down-cutting,
regular size,
19cm

33896.19

curettes - dissectors

Instruments @ 50% · Inserts to Scale



elevators

Instruments @ 50% · Inserts to Scale




FREER
Separator and Elevator
semi-sharp,
19.5cm

33461.19



FREER
Septum Elevator
strongly curved,
sharp/blunt tips,
19.5cm

33464.19



Ball-end Elevator and
Plugger, double-ended,
19.5cm

33477.19



GORNEY
Suction Elevator
obturator, blunt tip,
standard pattern,
20cm

33861.20




GORNEY
Suction Elevator
obturator, square tip,
semi-sharp,
20cm

33862.20


elevators

Instruments @ 50% · Inserts to Scale



DRAF
Micro Elevators
sharp, round handle, slightly
curved, 25cm

- A 33914.01** 1mm wide
- B 33914.03** 3mm wide







DRAF
Micro Elevators
sharp, round handle, 25cm

- A 33915.03** Strongly curved,
3mm wide
- B 33916.04** Angled,
4mm wide

cannulas - irrigators

Instruments @ 50% · Inserts to Scale





Goldenace P/L ACN 070 452 943
ABN 28 318 606 161

Head Office
116 Cutler Road, Jandakot WA 6164

Ph 08 9417 9545
Freecall 1800 756 757

quotes@inkasurgical.com.au
www.inkasurgical.com.au