

ENT

Paranasal Sinus Surgery



DISTRIBUTED BY:



SURGICAL ASIA Sdn Bhd (1272656-H)

Suppliers of Fine Quality Instruments

No. 8, Lebuhr Perusahaan Klebang 9A,
Kawasan Perusahaan IGB, 31200 Chemor, Perak, Malaysia

Phone: +605 210 8080 • www.inkasurgical.asia



Exclusive Distributor

OBEX Medical Ltd

0800 656 239

www.inkasurgical.co.nz



we thank you for your support

Disclaimer: Product availability is subject to regulatory approval in individual markets. Not all products are available in all markets. Contact us to find out more.

BRC-00155

contents

Alligator Forceps - Scissors

Cup Forceps


Nasal Cutting - Punch Forceps

Hook - Round - Sickle Knives

Curettes - Elevators

PARANASAL SINUS SURGERY

Delicate Scissors/Forceps for Paranasal Sinuses and Anterior Skull Base



BELLUCCI


Long Micro Alligator Scissors
straight blades 5mm
delicate

31650.13/13 shaft 13cm (5")

31650.13/15 shaft 15cm (6")

31650.13/19 shaft 19cm (7½")

also known as DRAF Ethmoid Scissors



HARTMANN

Long Alligator Forceps
serrated


30420.14 shaft 14cm (5½")

30420.16 shaft 16cm (6¼")

also known as DRAF Ethmoid Forceps



PARANASAL SINUS SURGERY

Delicate Sinus Forceps for Paranasal Sinuses and Anterior Skull Base



HARTMANN
Round Cup Forceps
ø 2.5mm
straight

30445.25/12 long shaft
12.5cm (5")
30445.25/15 X-long shaft
15cm (6")




WULLSTEIN
Micro Round Cup Forceps
ø 1.0mm
straight


31588.10/13 long shaft
13cm (5")
31588.10/15 X-long shaft
15cm (6")

PARANASAL SINUS SURGERY

Delicate Sinus Forceps for Paranasal Sinuses and Anterior Skull Base




STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm
long shaft 13cm (5")
33257.25/13 straight
33257.60/13 60° up-curved
33257.90/13 90° up-turned



STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm
X-long shaft 16cm (6¼")
33257.25/16 straight
33257.60/16 60° up-curved
33257.90/16 90° up-turned


also known as
DRAF Ethmoid Cup Forceps



STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm
XX-long shaft 19cm (7½")
33257.25/19 straight
33257.60/19 60° up-curved
33257.90/19 90° up-turned

PARANASAL SINUS SURGERY


Sinus Cup Forceps - Round Cup Forceps



Sinus Cup Forceps
with round cups
ø 3.5mm
shaft 11.5cm (4½")


33241.35

also known as HILDYARD Post-Nasal Biopsy Forceps




PARANASAL SINUS SURGERY

Nasal Cutting Forceps - Nasal Punch Forceps




WATSON WILLIAMS
Nasal Forceps
through-cutting
3mm oval cut
shaft 14cm (5½")
33301.03 straight
33302.03 up-curved

also known as WIGAND



MYLES-WIGAND
Nasal Punch Forceps
with flat jaws, serrated
through-cutting
2.5 x 11.5mm cut
shaft 14cm (5½")
33307.25 straight
33308.25 30° up-curved


also known as HERBERHOLD



GRUENWALD
Nasal Cutting Forceps
with Luer-Lock cleaning channel
through-cutting
working length 13cm (5")
straight
33338.20 2.0mm cut
33338.25 2.5mm cut
33338.30 3.0mm cut
up-curved
33339.20 2.0mm cut
33339.25 2.5mm cut
33339.30 3.0mm cut

PARANASAL SINUS SURGERY

Sickle Knives - Round Knives - Micro Elevators - Micro Curettes



original
size




DRAF
Sickle Knives
pointed tip
19cm (7½")

33901.19 straight shaft
33906.19 bayonet shaft



DRAF
Flap Knives
oval blade
19cm (7½")


33902.19 straight shaft
33907.19 bayonet shaft



DRAF
Round Knives
45° angled
19cm (7½")

straight shaft
33903.15 ø 1.5mm
33903.20 2.0mm

bayonet shaft
33908.15 ø 1.5mm
33908.20 2.0mm



33914.01

33914.03

33915.03



33916.04



DRAF
Micro Elevators
sharp
with round handle
25cm (10")

33914.01 slightly curved
1mm wide

33914.03 slightly curved
3mm wide

33915.03 strongly curved
3mm wide

33916.04 angled
4mm wide



Hook Knife
retrograde
straight
16.5cm (6½")
33909.16



HOUSE
Micro Curette
with oval cups
size 2
17cm (6¾")
31210.02



Goldenace P/L ACN 070 452 943
ABN 28 318 606 161

Head Office
116 Cutler Road, Jandakot WA 6164

Ph 08 9417 9545
Freecall 1800 756 757

quotes@inkasurgical.com.au
www.inkasurgical.com.au