

ENT

Rhinology



I N K A™

SURGICAL INSTRUMENTS

DISTRIBUTED BY:



SURGICAL ASIA Sdn Bhd (1272656-H)

Suppliers of Fine Quality Instruments

No. 8, Lebuhr Perusahaan Klebang 9A,
Kawasan Perusahaan IGB, 31200 Chemor, Perak, Malaysia

Phone: +605 210 8080 • www.inkasurgical.asia



Exclusive Distributor

OBEX Medical Ltd

0800 656 239

www.inkasurgical.co.nz



we thank you for your support

Disclaimer: Product availability is subject to regulatory approval in individual markets.
Not all products are available in all markets. Contact us to find out more.

BRC-00086

contents

Nasal Specula

Applicators - Probes

Septum Scissors

Snares

Dressing - Packing - Polypus Forceps

Nasal Cutting Forceps - Punches

Antrum Punches

Knives

Dissectors - Elevators - Probes - Raspatories

Chisels - Gouges - Osteotomes

Mallets

Rasps

Curettes

Retractors


Trocars - Cannulas

Suction Tubes

Needleholders


speculum

Instruments @ 50% unless stated · Inserts to Scale



THUDICUM
Nasal Specula
7cm (2¾)

- 33101.01** size 1 (11mm)
- 33101.02** size 2 (12mm)
- 33101.03** size 3 (15mm)
- 33101.04** size 4 (17mm)
- 33101.05** size 5 (19mm)
- 33101.06** size 6 (21mm)
- 33101.07** size 7 (23mm)







THUDICUM
Nasal Specula
insulated
7cm (2¾)

- 33104.01** size 1
- 33104.02** size 2
- 33104.03** size 3
- 33104.04** size 4
- 33104.05** size 5
- 33104.06** size 6
- 33104.07** size 7



for Diathermy



ST. CLAIR THOMSON
Nasal Specula
7.5cm (3")


- 33108.01** size 1 (38mm)
- 33108.02** size 2 (51mm)
- 33108.03** size 3 (64mm)
- 33108.04** size 4 (76mm)

speculum

Instruments @ 50% unless stated · Inserts to Scale

TIECK-HALLE
Nasal Specula
for infants
13.5cm (5¼")



33112.01 small
33112.02 large



VIENNA
Nasal Specula
14cm (5½")

33116.01 small
33116.02 medium
33116.03 large

standard
pattern



HARTMANN-HALLE
Nasal Specula
15.5cm (6")

33120.01 small
33120.02 medium
33120.03 large



BECKMANN
Nasal Speculum
curved to side
14.5cm (5¾")

33122.02




BECKMANN
Nasal Speculum
for infants
curved to side
13cm (5")

33122.13

speculum


Instruments @ 50% unless stated · Inserts to Scale



standard
pattern


COTTLE
Septum Specula
with thin tapered blades
14.5cm (5¾)

33125.01 35mm
33125.02 50mm
33125.03 75mm
33125.04 90mm





COTTLE
Septum Specula
with thin tapered blades
black coated
14.5cm (5¾)

33125.01 BL 35mm
33125.02 BL 50mm
33125.03 BL 75mm
33125.04 BL 90mm



COTTLE
Septum Specula
with lateral set screw
thin tapered blades
14.5cm (5¾)

33128.02 50mm
33128.03 75mm
33128.04 90mm





COTTLE
Septum Specula
with lateral set screw
thin tapered blades
black coated
14.5cm (5¾)

33128.02 BL 50mm
33128.03 BL 75mm
33128.04 BL 90mm


speculum

Instruments @ 50% unless stated · Inserts to Scale





speculum

Instruments @ 50% unless stated · Inserts to Scale



KILLIAN
Septum Specula
15cm (6")



33134.35 35mm
33134.50 50mm
33134.64 64mm
33134.75 75mm
33134.90 90mm



KILLIAN
Septum Specula
with **insulated** blades
and lateral set screw
15cm (6")

33137.50 50mm
33137.64 64mm
33137.75 75mm
33137.90 90mm

for diathermy



take apart
for easy cleaning

KILLIAN
Septum Specula
with **fiber optic light carrier**
and take-apart joint
15cm (6")

33136.40 40mm
33136.55 55mm
33136.75 75mm
33136.90 90mm
33136.98 Spare Fiber Optic Light Carrier
only


Triple Connector with adaptors
to fit light cables



ACMI
Brit.STANDARD




WOLF



STORZ
OLYMPUS

speculum



Instruments @ 50% unless stated · Inserts to Scale



KILLIAN-STRUYCKEN
Nasal Specula
with take-apart joint
15cm (6")

33142.01 left blade 75mm
right blade 45mm

33142.02 left blade 45mm
right blade 75mm




KILLIAN-STRUYCKEN
Septum Specula
slender blades
with take-apart joint
15cm (6")

33143.40 40mm

33143.55 55mm


33143.75 75mm

33143.90 90mm



probes - applicators

Instruments @ 50% unless stated · Inserts to Scale



JOBSON HORNE
Probe with ring
fine pattern
with knurled handle
and curved ring
18cm (7")
15659.18



JOBSON HORNE
Probe with ring
UK pattern
smooth round handle
with curved ring
18cm (7")
15658.18



FARRELL
Cotton Applicators
with round serrated tip
15667.14 14cm (5½")
15667.16 16cm (6¼")
15667.18 18cm (7")


also known as **LATHBURY**




FARRELL
Cotton Applicators
triangular serrated tip
15668.14 14cm (5½")
15668.16 16cm (6¼")
15668.18 18cm (7")

septum scissors


Instruments @ 50% unless stated · Inserts to Scale



CAPLAN
Septum Scissors
serrated blades
with angled shanks
21cm (8¼)
33386.21



CAPLAN
Septum Scissors
serrated blades
straight
21cm (8¼)
33386.22



CAPLAN-TC
Tungsten Carbide tipped
Septum Scissors
serrated blades
with angled shanks
21cm (8¼)
33387.21

gold plated handles

scissors


Instruments @ 50% unless stated · Inserts to Scale



HEYMANN
Nasal Scissors
with serrated blade
18cm (7")

33180.18

also known as LAKESIDE



HEYMANN-SC
SuperCut
Nasal Scissors
18cm (7")

33180.18 SC


black handles



HEYMANN-TC
Tungsten Carbide tipped
Nasal Scissors
with serrated blade
18cm (7")

33181.18

gold plated handles



HEYMANN
Turbinectomy Scissors
heavy pattern
with serrated blade
19.5cm (7 3/4)

33183.19



HEYMANN-SC
SuperCut
Turbinectomy Scissors
heavy pattern
19.5cm (7 3/4)

33183.19 SC

black handles




HEYMANN-TC
Tungsten Carbide tipped
Turbinectomy Scissors
heavy pattern
with serrated blade
19.5cm (7 3/4)

33184.19

gold plated handles


septum - turbinate scissors

Instruments @ 50% unless stated · Inserts to Scale



BECKER
Septum Scissors
serrated blades
with angled shanks
19cm (7½)
33384.19

GORNEY
Turbinate Scissors
long narrow blades
serrated
with angled shanks
21cm (8¼)
33390.21






GORNEY-TC
Tungsten Carbide tipped
Turbinate Scissors
long narrow blades
serrated
with angled shanks
21cm (8¼)
33391.21

gold plated handles


scissors

Instruments @ 50% unless stated · Inserts to Scale




snares

Instruments @ 50% unless stated · Inserts to Scale



KRAUSE
Nasal Snare
25cm (10")
33203.25





LANGE-WILDE
Nasal Snare
24cm (9½")
33205.24



Snare Wire
coil of 10m
ø 0.3mm (30G)
30620.03




Snare Wires, precut
10cm (4")
ø 0.3mm (30G)
33211.03 pack of 10




UNIVERSAL-TC
Tungsten Carbide tipped
Wire Cutting Scissor
one blade serrated
12cm (4¾")
10618.12
gold plated handles

forceps


Instruments @ 50% unless stated · Inserts to Scale




LUCAE
Nasal Dressing Forceps
14cm (5½")
12395.14 serrated
12396.14 1x2 teeth




JANSEN
Nasal Dressing Forceps
16cm (6¼")
12395.16 serrated
12396.16 1x2 teeth
also known as GRUENWALD




GRUENWALD
Nasal Dressing Forceps
19cm (7½")
12395.19 serrated
12396.19 1x2 teeth
also known as JANSEN




JANSEN-GRUENWALD
Nasal Dressing Forceps
serrated
21cm (8¼")
12395.21




WILDE
Nasal Dressing Forceps
angled
15cm (6")
33155.15
also known as TROELTSCH



TROELTSCH
Nasal Dressing Forceps
angled
17.5cm (7")
33155.17




TROELTSCH
Nasal Tissue Forceps
1x2 teeth
angled
15cm (6")
33156.15




TROELTSCH
Nasal Dressing Forceps
paediatric pattern
angled
11cm (4½")
30395.11
also known as WILDE

forceps


Instruments @ 50% unless stated · Inserts to Scale




HARTMANN
Nasal Dressing Forceps
18cm (7")
33160.18




HARTMANN
Nasal Dressing Forceps
20cm (8")
33163.20



WESTMACOTT
Nasal Dressing Forceps
delicate
20cm (8")
33167.20
also known as WEST




LISTER
Sinus Forceps
straight
13620.12 12.5cm (5")
13620.15 15.5cm (6")



TILLEY
Nasal Dressing Forceps
17.5cm (7")
13624.17

forceps


Instruments @ 50% unless stated · Inserts to Scale



NOYES
Nasal Polypus Forceps
shaft 11.5cm (4½")
33171.11 serrated jaws
33172.11 1x2 teeth



forceps

Instruments @ 50% unless stated · Inserts to Scale




standard
pattern


LUC
Nasal Turbinate Forceps
with mortised shanks
oval shaped
19cm (7½)
33223.01 small
33223.02 medium
33223.03 large




LUC
Nasal Turbinate Forceps
with mortised shanks
heart shaped
19cm (7½)
33226.01 small
33226.02 medium
33226.03 large




mortised shanks offer greater stability - no overlapping jaws



LUC
Ethmoid Forceps
with fenestrated jaws
light pattern
21cm (8¼)
33216.01 small (9mm)
33216.02 large (12mm)




IRWIN MOORE
Nasal Turbinate Forceps
with mortised shanks
20cm (8")
33219.20




BRUENINGS
Septum Forceps
with mortised shanks
19cm (7½)
33230.01 size 1 - 7mm
33230.02 size 2 - 9mm

forceps

Instruments @ 50% unless stated · Inserts to Scale




KNIGHT
Septum Forceps
19cm (7½")
33232.19




CRAIG
Septum Bone Cutting Forceps
18cm (7")
33234.18


also known as **VOMER Forceps**



CRAIG-DOMNICK
Septum Bone Cutting Forceps
curved to left
18cm (7")
33235.18




LEWIS
Septum Forceps
with grooved jaws
11mm ø
opening to left side
shaft 11.5cm (4½")
33237.11




KRESSNER
Turbinate Crushing Forceps
21cm (8¼")
33239.03 3mm wide
33239.04 4mm wide

septum forceps

Instruments @ 50% unless stated · Inserts to Scale




JANSEN-MIDDLETON
Septum Forceps
with compound action
spoon-shaped 4.8x14mm
18cm (7")
33376.18



JANSEN-STRUYCKEN
Septum Punch Forceps
with compound action
through-cutting
19cm (7½)
33380.03 small pattern
3.5mm cut
33380.04 regular pattern
4.5mm cut

forceps

Instruments @ 50% unless stated · Inserts to Scale



WEIL-BLAKESLEY
Sinus Forceps
fenestrated jaws
shaft 11.5cm (4½)

straight

33244.25 2.5mm bite

33244.30 3.0mm bite

33244.36 3.6mm bite

33244.40 4.0mm bite

33244.50 5.0mm bite

45° up-curved

33245.25 2.5mm bite

33245.30 3.0mm bite

33245.36 3.6mm bite

33245.40 4.0mm bite

33245.50 5.0mm bite


90° up-turned

33246.36 3.6mm bite

33246.40 4.0mm bite

33246.50 5.0mm bite


also known as BLAKESLEY-WILDE



BLAKESLEY
Nasal Forceps
with jaws 30° curved
3.6mm bite
shaft 11.5cm (4½)

33249.36 curved left

33250.36 curved right




also known as HILDYARD Post-Nasal Biopsy Forceps

Sinus Cup Forceps
with round cups
ø 3.5mm
shaft 11.5cm (4½)

33241.35


forceps

Instruments @ 50% unless stated · Inserts to Scale




WILDE
Ethmoid Forceps
fenestrated jaws
shaft 14cm (5½")
straight
21609.03 3mm bite
21609.04 4mm bite
21609.05 5mm bite
21609.06 6mm bite
up-curved
21612.03 3mm bite
21612.04 4mm bite
21612.05 5mm bite
21612.06 6mm bite

also known as WEIL-BLAKESLEY long pattern




TAKAHASHI
Ethmoid Forceps
shaft 12cm (4¾")
straight
21617.02 2.5mm bite
21617.04 4.0mm bite
up-curved
21618.02 2.5mm bite
21618.04 4.0mm bite



ST. BARTS
Ethmoid Forceps
3mm bite
shaft 11.5cm (4½")
21628.11 straight
21629.11 up-curved
21630.11 down-curved
also known as
ST. BARTHOLOMEWS
or BEAD Fore ps

forceps


Instruments @ 50% unless stated · Inserts to Scale



WIGAND
Nasal Forceps
60° up-curved
cup jaws 14x6mm
with rounded tip
shaft 11.5cm (4½")


33269.60

also known as
BLAKESLEY-WIGAND



HARTMANN
Round Cup Forceps
ø 2.5mm
straight

30445.25/12 long shaft
12.5cm (5")
30445.25/15 X-long shaft
15cm (6")







STRUEMPEL
Oval Cupped Forceps
fenestrated jaws 2.5mm
shaft 10cm (4")

33257.25 straight
33257.45 45° up-curved
33257.90 90° up-turned

also known as **STRUEMPEL-VOSS**


forceps

Instruments @ 50% unless stated · Inserts to Scale




forceps

Instruments @ 50% unless stated · Inserts to Scale



WEIL
Nasal Forceps
with oval cup jaws
5mm bite
shaft 11.5cm (4½")
33272.05 straight
33273.05 up-curved







GRUENWALD
Nasal Forceps
straight
shaft 12.5cm (5")
21672.02 2mm bite
21672.03 3mm bite
21672.04 4mm bite

The name changes according to the bite size. However, GRUENWALD has been the most commonly used. We use that name to make it less confusing. If you are not sure please talk to us


forceps

Instruments @ 50% unless stated · Inserts to Scale




forceps


Instruments @ 50% unless stated · Inserts to Scale



LUSK style
Mini Nasal Punch Forceps
with tubular shaft
through-cutting
1.5mm cut
13cm (5") working length
33318.15 straight
33318.45 up-curved




STRUYCKEN
Turbinate Forceps
through-cutting
1mm cut
shaft 11.5cm (4½)
33320.10



STRUYCKEN mod.
Turbinate Forceps
with improved jaws
through-cutting
1mm cut
shaft 11.5cm (4½)
33321.10

forceps - conchotomes


Instruments @ 50% unless stated · Inserts to Scale



MYLES-WIGAND
Nasal Punch Forceps
with flat jaws, serrated
through-cutting
2.5 x 11.5mm cut
shaft 14cm (5½)


33307.25 straight
33308.25 30° up-curved

also known as HERBERHOLD



HARTMANN
Conchotomes
with round jaws
through-cutting
shaft 11cm (4½)

33304.05 ø 5mm cut
33304.07 7mm cut
33304.09 9mm cut






SCHMEDEN
Nasal Cutting Forceps
with triangular jaws
through-cutting
shaft 11cm (4¼)

33315.06 6mm cut
33315.08 8mm cut
33315.10 10mm cut




suction forceps

Instruments @ 50% unless stated · Inserts to Scale



suction forceps


Instruments @ 50% unless stated · Inserts to Scale



through cutting forceps

Instruments @ 50% unless stated · Inserts to Scale

	outside width	-	inside cut
	3.3mm	-	2.0mm
	3.9mm	-	2.5mm
	4.6mm	-	3.0mm
	5.1mm	-	3.5mm
	5.7mm	-	4.0mm




WEIL-BLAKESLEY
Through cutting Nasal Forceps
with serrated lower jaw
shaft 11.5cm (4½)

straight
33261.20 2.0mm cut
33261.25 2.5mm cut
33261.30 3.0mm cut
33261.35 3.5mm cut
33261.40 4.0mm cut


up-curved
33264.20 2.0mm cut
33264.25 2.5mm cut
33264.30 3.0mm cut
33264.35 3.5mm cut
33264.40 4.0mm cut

also known as **BLAKESLEY**



WEIL-BLAKESLEY
Through cutting Nasal Forceps
with serrated lower jaw
curved jaws 3.6mm wide
2.5mm cut
shaft 11.5cm (4½)

33262.36 left curved
33265.36 right curved




WEIL-BLAKESLEY
Through cutting Nasal Forceps
with suction
shaft 11.5cm (4½)

straight
33263.25 2.5mm wide
33263.36 3.6mm wide
33263.48 4.8mm wide

up-curved
33266.25 2.5mm wide
33266.36 3.6mm wide
33266.48 4.8mm wide

nasal cutting forceps

Instruments @ 50% unless stated · Inserts to Scale




GRUENWALD
Nasal Forceps
through-cutting
shaft 11cm (4½)

33291.01 size 1 - 3.0mm cut


33291.02 size 2 - 3.5mm cut

33291.03 size 3 - 4.0mm cut



GRUENWALD mod.
Nasal Forceps
through-cutting
3.5mm wide
shaft 11cm (4½)

33292.35



MYLES
Nasal Cutting Forceps
through-cutting
shaft 12cm (4¾)

33294.01 size 1 - 2.0mm cut


33294.02 size 2 - 2.6mm cut

33294.03 size 3 - 3.0mm cut

33294.04 size 4 - 3.5mm cut

nasal cutting forceps

Instruments @ 50% unless stated · Inserts to Scale




WATSON WILLIAMS
Nasal Forceps

through-cutting
3mm oval cut
shaft 14cm (5½)

33301.03 straight

33302.03 up-curved

also known as WIGAND





MACKAY "Chomper"
Nasal Forceps

through-cutting
5.5mm wide
shaft 11.5cm (4½)

33305.55


nasal cutting forceps

Instruments @ 50% unless stated · Inserts to Scale




terrier punches

Instruments @ 50% unless stated · Inserts to Scale



OSTROM "Terrier"
Antrum Punches
backward cutting
with large oval cut 3.8mm
shaft 11cm (4½)

33342.38 down-curved
33343.38 up-curved




OSTROM "Terrier"
Antrum Punches
paediatric size
backward cutting
with small oval cut 2.0mm
shaft 11cm (4½)


33342.20 down-curved
33343.20 up-curved

backbiter punches


Instruments @ 50% unless stated · Inserts to Scale



OSTROM "Backbiter"
Antrum Punches
paediatric size
1.5mm cut
10cm (4") working length
33345.15 straight
33347.15 left-cutting
33348.15 right-cutting




OSTROM "Backbiter"
Antrum Punches
standard size
2.0mm cut
10cm (4") working length
33345.20 straight
33347.20 left-cutting
33348.20 right-cutting




OSTROM "Backbiter"
Antrum Punches
large size
2.5mm cut
10cm (4") working length
33345.25 straight
33347.25 left-cutting
33348.25 right-cutting

backbiter punches

Instruments @ 50% unless stated · Inserts to Scale






OSTROM "Backbiter"
Antrum Punches
upward cutting
10cm (4") working length
33350.15 1.5mm cut
33350.20 2.0mm cut
33350.25 2.5mm cut



OSTROM-STAMMBERGER
"Backbiter" Antrum Punch
backward cutting
2.5mm cut
shaft 11.5cm (4½")
33351.25
also known as American pattern

rotatable backbiter punches


Instruments @ 50% unless stated · Inserts to Scale




also known as TOBEY Micro Backbiter

ostrom punches


Instruments @ 50% unless stated · Inserts to Scale



BACKBITER
OSTROM Antrum Punches
paediatric size
with Luer-Lock cleaning channel
cut 1.5x5mm
10cm (4") working length
33359.15 up-cutting
33360.15 left-cutting
33361.15 right-cutting





BACKBITER
OSTROM Antrum Punches
with Luer-Lock cleaning channel
cut 2.5x7mm
10cm (4") working length
33359.25 up-cutting
33360.25 left-cutting
33361.25 right-cutting



OSTROM
Antrum Suction Punches
backward cutting
10cm (4") working length
33364.20 cut 2.0mm
shaft ø 4mm
33364.25 cut 2.5mm
shaft ø 5mm


sidebiting punches

Instruments @ 50% unless stated · Inserts to Scale




punches

Instruments @ 50% unless stated · Inserts to Scale




KERRISON
Antrum Punches
shaft 9cm (3½")

21762.35	3.5mm bite
21762.40	4.0mm bite
21762.45	4.5mm bite
21762.50	5.0mm bite
21762.60	6.0mm bite




KERRISON
Antrum and Sphenoid Punches
shaft 15cm (6")

40° forward cutting		90° up-cutting
21793.02	2mm	21795.02
21793.03	3mm	21795.03
21793.04	4mm bite	21795.04



CITELLI
Antrum Punches
through-cutting
shaft 7.5cm (3")

30570.20	2mm cut
30570.30	3mm cut






SEYMOUR JONES
Antrum Punch
4mm bite
15cm (6")

33394.15




punches

Instruments @ 50% unless stated · Inserts to Scale




punches

Instruments @ 50% unless stated · Inserts to Scale



kerrison punches

Instruments @ 50% unless stated · Inserts to Scale



KERRISON
Semi-detachable
Sphenoid Punches
shaft 18cm (7")

90° up-cutting

21910.01 1mm bite

21910.02 2mm bite

21910.03 3mm bite

21910.04 4mm bite


40° forward cutting

21921.01 1mm bite

21921.02 2mm bite

21921.03 3mm bite


21921.04 4mm bite




*easy to clean:
4 simple steps
to dismantle & re-assemble
quick & easy*

sphenoidal punches

Instruments @ 50% unless stated · Inserts to Scale





HAJEK-KOFLER
Sphenoidal Punch
360° rotatable
with fixation nut
3mm bite
shaft 18cm (7")
33397.03



MC KENTY
Sphenoidal Punch
through-cutting
2mm cut
360° rotatable
with fixation nut
shaft 15cm (6")
33400.02

knives

Instruments @ 50% unless stated · Inserts to Scale




dissectors


Instruments @ 50% unless stated · Inserts to Scale




ST.CLAIR THOMSON
Dissectors
16cm (6¼")
33451.16 sharp
33452.16 blunt




FREER
Dissectors
19cm (7½")
33456.19 sharp
33457.19 blunt




FREER
Dissector
sharp/blunt
18cm (7")
15643.18



FREER
Dissectors
long pattern
sharp/blunt
20cm (8")
15644.04 4.5mm wide
15644.06 6.5mm wide



FREER-USA
Dissectors
with square handle
sharp/blunt
19.5cm (7¾")
15646.04 4.5mm wide
15646.06 6.0mm wide




Mc DONALD
Dissector
blunt
19cm (7½")
15632.19

America n pattern



elevators

Instruments @ 50% unless stated · Inserts to Scale




dissectors - elevators

Instruments @ 50% unless stated · Inserts to Scale




dissectors - elevators - raspatory


Instruments @ 50% unless stated · Inserts to Scale



FREER
Separator and Elevator
semi-sharp
19.5cm (7¾")
33461.19



FREER
Nasal Raspatory
and Septum Knife
sharp
18cm (7")
33462.18



FREER-JOSEPH
Septum Elevator
and Raspatory
sharp
20cm (8")
33459.20



FREER
Septum Elevator
strongly curved
sharp/blunt
19.5cm (7¾")
33464.19




33467.19 33467.22

BALLENGER-HAJEK
Dissectors
sharp/blunt
33467.19 19.5cm (7¾")
33467.22 22.0cm (8¾")


elevators - raspatories

Instruments @ 50% unless stated · Inserts to Scale



5mm

JOSEPH
Periosteal Elevator
sharp
16cm (6 $\frac{3}{4}$)
20323.16



5mm

LEMPERT-JOSEPH
Periosteal Elevator
sharp
17.5cm (7")
20323.17




5mm

CONVERSE
Periosteal Elevator
sharp
18cm (7")
20327.18




6mm

Periosteal Elevator
with phenolic handle
sharp
18cm (7")
20387.18



8mm

COTTLE
Raspatory
sharp
19.5cm (7 $\frac{3}{4}$)
33483.19




9mm

JANSEN
Periosteal Elevator
sharp
19.5cm (7 $\frac{3}{4}$)
30524.09



6mm

ERLANGEN
Periosteal Elevator
sharp
17.5cm (7")
33484.17




6mm


JOSEPH
Nasal Elevator
sharp
18cm (7")
33487.18

elevators - osteotomes - gouges


Instruments @ 50% unless stated · Inserts to Scale




Septum Elevator
with angled blade
sharp
16.5cm (6½)
33489.16




WULLSTEIN-STORZ
Elevator
sharp
17cm (6¾)
31222.17




SEWALL
Antrum Elevator
semi-sharp
19.5cm (7¾)
33491.19




FREER
Septum Elevator
with square edge
sharp
17cm (6¾)
33492.17




WEST
Septum Osteotome
straight
16cm (6¼)
33496.16



WEST
Septum Osteotome
straight
19cm (7½)
33496.19




WEST
Septum Gouge
straight
16cm (6¼)
33499.16



WEST
Septum Gouge
straight
19cm (7½)
33499.19


chisels - gouges

Instruments @ 50% unless stated · Inserts to Scale




33503.16 **33504.16** 4mm

FREER
Septum Chisels
16cm (6¼")
33503.16 straight
33504.16 curved




FREER
Septum Chisel
19cm (7½")
33509.19


also known as Periodontome




FREER
Septum Gouge
straight
15cm (6")
33513.15



BALLENGER
Septum Gouges
with V-notched tip
15cm (6")
33517.04 4mm wide
33517.06 6mm wide
33517.08 8mm wide




KILLIAN
Septum Gouge
bayonet shank
17cm (6¾")
33521.15



KILLIAN-CLAUS
Septum Gouge
with V-notched tip
17cm (6¾")
33525.15



also known as TILLEY



KILLIAN-CLAUS mod.
Septum Gouge
with V-shaped tip
16.5cm (6½")
33525.16


mallets

Instruments @ 50% unless stated · Inserts to Scale



rasps



Instruments @ 50% unless stated · Inserts to Scale




Glabella Rasp
Tungsten Carbide tipped
down-cutting
21cm (8¼")
33537.21



blade dusted
with
Tungsten Carbide
particles




PEET-AUFRICHT
Glabella Rasp
with TC-dusted blade
21cm (8¼")
33540.21



Push & Pull
action



AUFRICHT
Glabella Rasp
up & down-cutting
20cm (8")
33544.20



up-cutting
Push
action

down-cutting
Pull
action



AUFRICHT mod.
Glabella Rasps
with fine teeth
20.5cm (8")
33548.20 up-cutting
33549.20 down-cutting
also known as American pattern



up-cutting down-cutting

standard
pattern



AUFRICHT
Glabella Rasps
with coarse teeth
21cm (8¼")
33548.21 up-cutting
33549.21 down-cutting


rasps - otoabraders

Instruments @ 50% unless stated · Inserts to Scale




33553.18 **33554.18**

LEWIS
Nasal Rasps
down-cutting
18cm (7")
33553.18 fine teeth
33554.18 coarse teeth



33557.18 **33558.18**

LEWIS-TC
Tungsten Carbide tipped
Nasal Rasps
down-cutting
18cm (7")
33557.18 fine teeth
33558.18 coarse teeth
gold plated handles



DINGMAN
Otoabraders
15.5cm (6")
33532.15 curved left
33533.15 curved right
also known as Cartilage Abraders




WATSON WILLIAMS
Antrum Rasp
UK pattern
heavy
17cm (6¾")
33574.17




WATSON WILLIAMS
Antrum Rasp
German pattern
23cm (9")
33576.23

rasps - chisels - burs

Instruments @ 50% unless stated · Inserts to Scale




BRAWLEY
Antrum Rasp
20cm (8")
33579.20




33584.19 **33585.19**

standard
pattern

WIENER
Antrum Rasps
19.5cm (7¾")
33584.19 trocar tip
33585.19 blunt tip




WIENER
Antrum Rasp
with **small trocar tip**
19.5cm (7¾")
33584.20



FAULKNER
Antrum Trocar Chisel
9mm wide
15.5cm (6")
33588.15



TILLEY "Harpoon"
Antrum Trocar
17cm (6¾")
33598.17



6mm



7mm

8mm

TILLEY
Antrum Burs
17cm (6¾")
33601.06 ø 6mm
33601.07 7mm
33601.08 8mm



rasps - curettes

Instruments @ 50% unless stated · Inserts to Scale




curettes

Instruments @ 50% unless stated · Inserts to Scale





FAULKNER
Antrum Curette
double-ended
21.5cm (8½)
33643.21



COAKLEY
Antrum Curettes
17cm (6¾)

- 33652.01** size 1 - oval, light curve
- 33652.02** size 2 - oval, strong curve
- 33652.03** size 3 - oval, acute curve
- 33652.04** size 4 - small, light curve
- 33652.05** size 5 - small, strong curve
- 33652.06** size 6 - triangular, curved



MACKIE
Ethmoid Curettes
double-curved
19.5cm (7¾)


- 33648.06** 6mm wide
- 33648.08** 8mm wide
- 33648.09** 9mm wide

retractors

Instruments @ 50% unless stated · Inserts to Scale




LANGENBECK
Retractors have
ROUNDED blades




LANGENBECK
Retractors
21.5cm (8½")

- 15241.22** 20x7mm wide
- 15241.30** 30x10mm wide
- 15241.31** 30x14mm wide
- 15241.40** 40x10mm wide





LANGENBECK-KOCHER
Retractors have
FLAT blades




LANGENBECK-KOCHER
Retractors
22cm (8¾")


- 15245.25** 25x8mm wide
- 15245.35** 35x8mm wide
- 15245.36** 35x11mm wide
- 15245.37** 35x15mm wide
- 15245.40** 40x11mm wide
- 15245.55** 55x11mm wide



15198.14





15198.20




15198.26

MIDDELDORPF
Cheek Retractors


- 15198.14** 14x17mm wide
21.5cm (8½")
- 15198.20** 20x22mm wide
22.5cm (9")
- 15198.26** 26x30mm wide
23.5cm (9¼")




17mm



22mm





30mm



BRUENINGS-BRUNK
Antrum Retractor
toothed inside blade
45x15mm wide
19cm (7½")
33608.19

retractors

Instruments @ 50% unless stated · Inserts to Scale





retractors




Instruments @ 50% unless stated · Inserts to Scale




FISCH
Retractors
with self-arresting arms
3x3 prongs
17cm (6¾")
30722.17 semi-sharp
30723.17 blunt



30722.17
30723.17





PERKINS
Antrum Retractors
with 3 semi-sharp prongs
and solid grooved blade
with curved arms
13cm (5")
33615.13 solid blade left side
33616.13 solid blade right side



PERKINS
Antrum Retractors
with 3 semi-sharp prongs
and solid blade with pins
with curved arms
13cm (5")
33618.13 solid blade left side
33619.13 solid blade right side


trocars

Instruments @ 50% unless stated · Inserts to Scale




cannulas

Instruments @ 50% unless stated · Inserts to Scale




suction tubes - probes

Instruments @ 50% unless stated · Inserts to Scale



nasal suction tubes

Instruments @ 50% unless stated · Inserts to Scale




standard pattern

105mm

FRAZIER
Suction Tubes
with finger valve
and stilette
working length 105mm
18cm (7")


- 15705.03** FG 3 (1.0mm)
- 15705.04** FG 4 (1.2mm)
- 15705.05** FG 5 (1.7mm)
- 15705.06** FG 6 (2.0mm)
- 15705.07** FG 7 (2.3mm)
- 15705.08** FG 8 (2.7mm)
- 15705.09** FG 9 (3.0mm)
- 15705.10** FG 10 (3.3mm)
- 15705.12** FG 12 (4.0mm)
- 15705.15** FG 15 (5.0mm)
- 15705.99** Spare Stilettes only



100mm

FRAZIER
Suction Tubes
calibrated 1 - 9cm
with finger valve
and stilette
working length 100mm
18cm (7")

- 33700.06** FG 6 (2.0mm)
- 33700.07** FG 7 (2.3mm)
- 33700.08** FG 8 (2.7mm)
- 33700.09** FG 9 (3.0mm)
- 33700.10** FG 10 (3.3mm)
- 33700.12** FG 12 (4.0mm)
- 33700.14** FG 14 (4.7mm)
- 15705.99** Spare Stilettes only




100mm

FRAZIER
Suction Tubes
ball point tipped
with finger valve
and stilette
working length 100mm
18cm (7")

- 33702.08** FG 8
- 33702.10** FG 10
- 33702.12** FG 12
- 33702.15** FG 15
- 15705.99** Spare Stilettes only

when ordering please specify stilettes (obturators) please indicate FG size


Outer ø	FG	mm
	6	2.0
	7	2.3
	8	2.7
	9	3.0
	10	3.3
	12	4.0
	15	5.0



130mm

FRAZIER
Suction Tubes
with finger valve
and stilette
long pattern
working length 130mm
20cm (8")

- 15704.05** FG 5
- 15704.06** FG 6
- 15704.07** FG 7
- 15704.08** FG 8
- 15704.09** FG 9
- 15704.10** FG 10
- 15704.12** FG 12
- 15704.15** FG 15
- 15704.99** Spare Stilettes only



100mm





33709.10 **33710.10**

FRAZIER
Suction Tubes
with 45° curved tip
FG 10 (3.3mm)
with finger valve
and stilette
working length 100mm
18cm (7")

33709.10 curved to left
33710.10 curved to right
15705.99 Spare Stilettes only


nasal suction tubes

Instruments @ 50% unless stated · Inserts to Scale



nasal suction tubes


Instruments @ 50% unless stated · Inserts to Scale



ANDREWS-PYNCHON mod.
Small Suction Tube
with **removable tip**
23cm (9")

15733.23

15733.99 Spare Tip only




MAGILL
Suction Tubes
with stilette
chrome plated
21cm (8 3/4)

15701.01 Ø 3.0mm - small

15701.02 3.5mm - medium

15701.03 4.0mm - large

15701.99 Spare Stilettes
only




BABY-YANKAUER
Small Suction Tube
22cm (8 3/4)

15737.21

15737.99 Spare Tip only

nasal suction tubes - irrigators

Instruments @ 50% unless stated · Inserts to Scale



HOUSE
Suction Irrigators
with finger valve
and stilettes
working length 90mm
16cm (6 $\frac{1}{4}$)

15724.23 2.3/1.7mm (FG 7/5)
15724.25 2.5/2.0mm (FG 8/6)
15724.33 3.3/2.3mm (FG 10/7)
15724.40 4.0/2.7mm (FG 12/8)
suction/irrigation

FUKUSHIMA style
"Teardrop" Suction Tubes
tapered
malleable stainless steel
hub length 100mm
total length 16.5cm (6 $\frac{1}{2}$)


15727.04 FG 4 (1.3mm)
15727.05 FG 5 (1.7mm)
15727.06 FG 6 (2.0mm)
15727.07 FG 7 (2.3mm)
15727.08 FG 8 (2.7mm)
15727.09 FG 9 (3.0mm)
15727.10 FG 10 (3.3mm)
15727.12 FG 12 (4.0mm)

FUKUSHIMA style
"Teardrop" Suction Tubes
tapered
malleable stainless steel
hub length 115mm
total length 18cm (6 $\frac{1}{2}$)

15728.04 FG 4
15728.05 FG 5
15728.06 FG 6
15728.07 FG 7
15728.08 FG 8
15728.09 FG 9
15728.10 FG 10
15728.12 FG 12

needleholders


Instruments @ 50% unless stated · Inserts to Scale



Micro Needleholders
with **TC-dusted jaws**
very delicate
with black handles

14100.15 15.0cm (6")
14100.16 16.5cm (6½")


black handles



CRILE-WOOD-TC
Tungsten Carbide tipped
Needleholder
delicate
15cm (6")

14138.15

gold plated handles




16cm

18cm

CRILE-WOOD-TC
Tungsten Carbide tipped
Needleholders
delicate

14138.16 16cm (6¼")
14138.18 18cm (7")

gold plated handles




16cm

18cm

MAYO-HEGAR-TC
Tungsten Carbide tipped
Needleholders

14161.16 16cm (6¼")
14161.18 18cm (7")


gold plated handles



OLSEN-HEGAR-TC
Tungsten Carbide tipped
Needleholder with Scissors
17cm (6¾")

14201.17

gold plated handles



PARELL-TC
Intra-nasal Needleholder
Tungsten Carbide tipped
with recessed jaws
angled to side
18cm (7")

14155.18

gold plated handles



Goldenace P/L ACN 070 452 943
ABN 28 318 606 161

Head Office
116 Cutler Road, Jandakot WA 6164

Ph 08 9417 9545
Freecall 1800 756 757

quotes@inkasurgical.com.au
www.inkasurgical.com.au